



Recombinant Mouse CD40 ligand (Cd40lg)

Product Code	CSB-CF004937MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P27548
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MIETYSQPSPRSVATGLPASMKIFMYLLTVFLITQMIGSVLFAVYLHRRLDKVEE EVNLHEDFVFIKKLKRCNKGEGLSLLNCEEMRRQFEDLVKDITLNKEEKENS FEMQRGDEDPQIAAHVVSEANSNAASVLQWAKKGYYTMKSNLVMLENGKQL TVKREGLYYVYTQVTFCNREPSSQRPFIVGLWLKPSSGSERILLKAANTHSS SQLCEQQSVHLGGVFELQAGASVFVNVTEASQVIHRVGFSSFGLLKL
Source	in vitro E.coli expression system
Target Names	Cd40lg
Protein Names	Recommended name: CD40 ligand Short name= CD40-L Alternative name(s): T-cell antigen Gp39 TNF-related activation protein Short name= TRAP Tumor necrosis factor ligand superfamily member 5 CD_antigen= CD154 Cleaved into the follow
Expression Region	1-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.