



Recombinant Human Syntaxin-7 (STX7)

Product Code	CSB-CF022903HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15400
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	SYTPGVGGDPAQLAQRISSENIQKITQCSVEIQRTLNLGTPQDSPELRQQQLQQ KQQYTNQLAKETDKYIKEFGSLPTTPSEQRQRKIQKDRLVAEFTTSLTNFQKV QRQAAEREKEFVARVRASSRVSGSFPEDSSKERNLVSWEQSTQPQVQVQDE EITEDDLRLIHERESSIRQLEADIMDINEIFKDLGMMIHEQGDVIDSIEANVENAE VHVQQANQQLSRAADYQRKSRKTLCIILILVIGVAIISLIWGLNH
Source	in vitro E.coli expression system
Target Names	STX7
Protein Names	Recommended name: Syntaxin-7
Expression Region	2-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.