



Recombinant Human Tissue factor (F3)

Product Code	CSB-CF007928HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13726
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	SGTTNTVAAYNLTWKSTNFKTILEWEPKPVNQVYTVQISTKSGDWKSKCFYTT DTECDLTDEIVKDVKQTYLARVFSYPAGNVESTGSAGEPLYENSPEFTPYLET NLGQPTIQSFEQVGTKVNVTVEDERTLVRRNNTFLSLRDVFGKDLYTLYYWK SSSSGKKTAKTNTNEFLIDVDKGENYCFSVQAVIPSRTVNRKSTDSPVECMGQ EKGEFREIFYIIGAVV FVVIILVILLAISLHKCRKAGVGGQSWKENSPLNVS
Source	in vitro E.coli expression system
Target Names	F3
Protein Names	Recommended name: Tissue factor Short name= TF Alternative name(s): Coagulation factor III Thromboplastin CD_antigen= CD142
Expression Region	33-295
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.