



Recombinant Human Receptor-transporting protein 1 (RTP1)

Product Code	CSB-CF020577HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59025
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MRIFRPWRLRCPALHLP SLSVFS LRWKLPSLTTDETMCKSVTTDEWKKVFYEK MEEAKPADSWDLIIDPNLKHNVLS PGWVKQYLELHASGRFHCSWCWHTWQSP YVVILFHMFLDRAQRAGSVRMRVFKQLCYECGTARLDESSMLEENIEGLVDNL ITSLREQCYGERGGQYRIHVASRQDNRRHRGEFCEACQEGIVHWKPSEKLLLE EEATTYTFSRAPSPTKSQDQTGSGWNFCSIPWCLFWATVLLLIY LQFSFRSSV
Source	in vitro E.coli expression system
Target Names	RTP1
Protein Names	Recommended name: Receptor-transporting protein 1
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.