



Recombinant Rabbit Translocon-associated protein subunit alpha (SSR1)

Product Code	CSB-CF022717RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53815
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	RGGPGGSLAEAQDLSEDEETVEDSVIEDEDDEAEVEEDEPTDLAEDREEEDV SGEPEASPSADTTILFVKGEDFPANNIVRFLVGFTNKGTEDFIVESLDASFRYP QDYQFYIQNFTALPLNTIVPPQRQATFEYSFIPAEPMGGRPFGLVINLNYKDFN GNVFQDAVFNQTVTIEREDGLDGQTIFMYMSLAGLGLLVVGLHQLLESNRN KRPIQKVEMGTSSQNDVDMSWIPQETLNQINKASPRRLPRKRPQKRSSVGSDE
Source	in vitro E.coli expression system
Target Names	SSR1
Protein Names	Recommended name: Translocon-associated protein subunit alpha Short name= TRAP-alpha Alternative name(s): Signal sequence receptor subunit alpha Short name= SSR-alpha
Expression Region	22-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.