



Recombinant Rabbit T-lymphocyte activation antigen CD80 (CD80)

Product Code	CSB-CF004959RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42070
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	GISQVTKSVKEMAALSCDYNISIDELARMRIYWQKDQQMVLISIISGQVEVWPEY KNRTFPDIINLSLMILALRLSDKGYTCVVQKNENGSFRREHLTSVTLIRADF PVPSITDIGHPDPNVKRIRCSASGGFPEPRLAWMEDGEELNAVNTTVDQDLDT ELYSVSSELDNFVTNNHSIVCLIKYGELSVSQIFPWSKPKQEPPIDQLPFWVIIP VSGALVLTAVVLYCLACRHVARWKRTRRNEETVGTERLSPIYLGSAQSSG
Source	in vitro E.coli expression system
Target Names	CD80
Protein Names	Recommended name: T-lymphocyte activation antigen CD80 Alternative name(s): Activation B7-1 antigen CD_antigen= CD80
Expression Region	33-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.