



# Recombinant Mouse Tumor necrosis factor ligand superfamily member 6 (Faslg)

<b>Product Code</b>	CSB-CF008434MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P41047
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	MQQPMNYPQPQIFWVDSSATSSWAPPGSVFPCPSCGPRGPDQRRPPPPPP PVSPLPSPSQPLPLPPLTPLKDKDHNTNLWLPVVFMMVLVALVGMGLGMYQLF HLQKELAELEFTNQSLKLVSSFQKQIANPSTPSEKKEPRSVAHLTGNPHSR SIP LEWEDTYGTALISGVKYYKGGGLVINETGLYFVYSKVYFRGQSCNNQPLNHKVY MRNSKYPEDLVLMEEKRLNYCTTGQIWAHSSYLGAVFNLTSAHLYVNISQLS LINFEEKTFFGLYKL
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Faslg
<b>Protein Names</b>	Recommended name: Tumor necrosis factor ligand superfamily member 6 Alternative name(s): CD95 ligand Short name= CD95-L Fas antigen ligand Short name= Fas ligand Short name= FasL CD_antigen= CD178 Cleaved into the following 4 chains: 1
<b>Expression Region</b>	1-279
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.