



Recombinant Human Low affinity immunoglobulin gamma Fc region receptor II-a (FCGR2A)

Product Code	CSB-CF008540HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P12318
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	QAAAPPKAVLKLEPPWINVLQEDSVTLTCQGARSPESDSIQWFHNGNLIPTHT QPSYRFKANNNDSGEYTCQTGQTSLSDPVHLTVLSEWLVLQTPHLEFQEGETI MLRCHSWKDKPLVKVTFQNGKSQKFSHLDPTFSIPQANHSHSGDYHCTGNI GYTLFSSKPVTITVQVPSMGSSSPMGIIVAVVIATAVAIVAIVVALIYCRKKRIS ANSTDPVAAAQFEPGRQMIAIRKRQLEETNNDYETADGGYMTLNPRAPTDD DKNIYLTLPNDHVNSNN
Source	in vitro E.coli expression system
Target Names	FCGR2A
Protein Names	Recommended name: Low affinity immunoglobulin gamma Fc region receptor II-a Short name= IgG Fc receptor II-a Alternative name(s): CDw32 Fc-gamma RII-a Short name= Fc-gamma-RIIa Short name= FcRII-a CD_antigen= CD32
Expression Region	34-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.