



Recombinant Human ICOS ligand (ICOSLG)

Product Code	CSB-CF010958HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75144
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	DTQEKEVRAMVGSDELSCACPEGSRFDLNDVYVYWQTSESKTVVTYHIPQN SSELENVDSRYRNRALMSPAGMLRGDFSLRFLFNVTPQDEQKFHCLVLSQSLGF QEVLSVEVTLHVAANFSPVVSAPHSPSQDELFTFTCTTSINGYPRPNVYWINKT DNSLLDQALQNDTVFLNMRGLYDVVSVLRIARTPSVNIGCCIEVLLQQNLTVG SQTGNDIGERDKITENPVSTGEKNAATWSILAVLCLLVVAVAIGWVCRDRCLQ HSYAGAWAVSPETELTGHV
Source	in vitro E.coli expression system
Target Names	ICOSLG
Protein Names	Recommended name: ICOS ligand Alternative name(s): B7 homolog 2 Short name= B7-H2 B7-like protein GI50 B7-related protein 1 Short name= B7RP-1 CD_antigen= CD275
Expression Region	19-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.