



Recombinant Mouse T-lymphocyte activation antigen CD86 (Cd86)

Product Code	CSB-CF004965MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42082
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	VSVETQAYFNGTAYLPCPFTKAQNISSLVFWQDQQLVLYEHYLGTEKLD SVNAKYLGRTSFDRNNWTLRLHNVQIKDMGSYDCFIQKKPPTGSIILQQTLTEL SVIANFSEPEIKLAQNVTGNSGINLTCTSKQGHPKPKKMYFLITNSTNEYGDNM QISQDNVTELFISISNLSLSFPDGVWHMTVVCVLETESMKISSKPLNFTQEFPS PQTYWKEITASVTVALLLVMLLIIVCHKKPNQPSRPSNTASKLERDSNADRETIN LKELEPQIASAKPNAE
Source	in vitro E.coli expression system
Target Names	Cd86
Protein Names	Recommended name: T-lymphocyte activation antigen CD86 Alternative name(s): Activation B7-2 antigen Early T-cell costimulatory molecule 1 Short name= ETC-1 CD_antigen= CD86
Expression Region	24-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.