



Recombinant *Cricetulus griseus* Syndecan-1 (SDC1)

Product Code	CSB-CF020888DXU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47951
Product Type	Transmembrane Protein
Immunogen Species	<i>Cricetulus griseus</i> (Chinese hamster) (<i>Cricetulus barabensis griseus</i>)
Sequence	QIMAVNVPPEDQDGS GDDSDNFSGSGTGALPDITLSRQTSSTLKDVWLLTATP TAPEPTSRDTEATFTSILPAGEKPGEGEPVLIAEVDTSSTTW DKELEVTTTRPRE TTQLLVTHR VSTARATTAQAPVTSHPHRDVQPGLHETLAPTAPGQPDHQPPS GGTSVIKEVAEDGATNQLPTGEGSGEQDFTFETSGENTAVAAIEPDQRNQPP VDEGATGASQGLLDRKEVLGGVIAGGLVGLIFAVCLVGFMLYRMKKKDEGSYS LEEPKQANGGAYQKPTKQEEFYA
Source	in vitro E.coli expression system
Target Names	SDC1
Protein Names	Recommended name: Syndecan-1 Short name= SYND1 Alternative name(s): CD_antigen= CD138
Expression Region	23-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.