



Recombinant Bovine Syntaxin-1A (STX1A)

Product Code	CSB-CF022895BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32850
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MKDRTQELRTAKDSDDDDVTVTVDRDRFMDEFFEQVEEIRGFIDKISENVEE VKRKHSAILASPNPDEKTKEELEELMSDIKKTANKVRSKLSIEQSIEQEEGLNR SSADLRIRKTQHSTLSRKFVEVMSEYNATQSDYRERCKGRIQRQLEITGRITTT SEELEDMLESGNPAIFASGIIMDSSISKQALSEIETRHSEIIKLENSIRELHDMFM DMAMLVESQGEMIDRIEYNVEHSVDYVERAVSDTKKAVKYQSKARRKKIMIVIC CVVLGIVIASTFGGIFG
Source	in vitro E.coli expression system
Target Names	STX1A
Protein Names	Recommended name: Syntaxin-1A Alternative name(s): Neuron-specific antigen HPC-1 Synaptotagmin-associated 35 kDa protein Short name= P35A
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.