



Recombinant Human Syntaxin-2 (STX2)

Product Code	CSB-CF022898HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32856
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MRDRLPDLTACRKNDDGDTVVVVEKDHFMDDFFHQVEEIRNSIDKITQYVEEV KKNHSIILSAPNPEGKIKEELEDLNKEIKKTANKIRAKLKAIEQSFDQDESGNRTS VDLRIRRTQHSVLSRKFVEAMA EYNEAQLFRERSKGR IQRQLEITGR TTTDDE LEEMLESGKPSIFTSDIISDSQITRQALNEIESRHKDIMKLETSIRELHEMFMDMA MFVETQGEMINNIERNVMNATDYVEHAKEETKKAIKYQSKARRKKWIIIVSVV LVAIIALIIGLSVGK
Source	in vitro E.coli expression system
Target Names	STX2
Protein Names	Recommended name: Syntaxin-2 Alternative name(s): Epimorphin
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.