



Recombinant Mouse Syntaxin-1B (Stx1b)

Product Code	CSB-CF022896MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61264
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MKDRTQELRSAKDSDEEEVVHVD RDHFMDEF FEQVEEIRGCIEKLS EDVEQ VKKQHSAILAAPNPDEKTKQELEDLTADIKKTANKVRSKLKAIEQSIEQEEGLNR SSADLRIRKTQHSTLSRKFVEVMTEYNATQSKYRDRCKDRIQRQLEITGR TTTN EELEDMLESGKLAIFTDDIKMDSQMTKQALNEIETRHNEI IKLETSIRELHDMFV DMAMLVESQ GEMIDRIEYNVEHSVDYVERAVSDTKKAVKYQSKARRKKIMIIIC CVVLGVVLASSIGGTLGL
Source	in vitro E.coli expression system
Target Names	Stx1b
Protein Names	Recommended name: Syntaxin-1B
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.