



Recombinant Human Epithelial cell adhesion molecule (EPCAM)

Product Code	CSB-CF007717HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P16422
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	QEEVCENYKLAVNCFVNNNRQCQCTSVGAQNTVICSKLAAKCLVMKAEMN GSKLGRRAKPEGALQNNNDGLYDPDCDESGLFKAKQCNGTSMCWCVNTAGV RRTDKDTEITCSERVRTYWIIIELKHKAREKPYDSKSLRTALQKEITTRYQLDPK FITSILYENNVITIDLQNSSQKTQNDVDIADVAYYFEKDKVKGESLFHSHKMDLT VNGEQLDLDPGQTLIYYVDEKAPEFSMQGLKAGVIAVIVVVVIAVVAGIVVLVIS RKKRMAKYEKAEIKEMGEMHRELNA
Source	in vitro E.coli expression system
Target Names	EPCAM
Protein Names	Recommended name: Epithelial cell adhesion molecule Short name= Ep-CAM Alternative name(s): Adenocarcinoma-associated antigen Cell surface glycoprotein Trop-1 Epithelial cell surface antigen Epithelial glycoprotein Short name= EGP Epithe
Expression Region	24-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.