



Recombinant Mouse Asialoglycoprotein receptor 2 (Asgr2)

Product Code	CSB-CF002208MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P24721
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MEKDCQDIQQLDSEENDHQLSGDDEHGSHVQDPRIENPHWKGQPLSRPFPQ RLCSTFRLSLLALAFNILLVVICVSSQSIQLQEEFRTLKETFSNFSSSTLMEFG ALDTLGGSTNAILTSWLAQLEEKQQQLKADHSTLLFHLKHFPMDLRTLTCQLAY FQSNNGTECCPVNWVEFGGSCYWFSRDGLTWAEADQYCQLENAHLLVINSRE EQDFVVKHRSQFHIWIGLTDRDGSWKWVDGTDYRSNYRNWAFTQPDNWQG HEQGGGEDCAEILSDGHWNDNFCQQVNRWVCEKRRNITH
Source	in vitro E.coli expression system
Target Names	Asgr2
Protein Names	Recommended name: Asialoglycoprotein receptor 2 Short name= ASGP-R 2 Short name= ASGPR 2 Alternative name(s): Hepatic lectin 2 Short name= HL-2 Short name= mHL-2
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.