



Recombinant Pan troglodytes Sodium/potassium-transporting ATPase subunit beta-1 (ATP1B1)

Product Code	CSB-CF002326EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5A6J8
Product Type	Transmembrane Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Sequence	MARGKAKEEGSWKKFIWNSEKKEFLGRTGGSWFKILLFYVIFYGCLAGIFIGTI QVMLLTISEFKPTYQDRVAPPGLTQIPQIQKTEISFRPNDPKSYEAYVLNIVRFL EKYKDSAQRDDMIFEDCGDVPSEPKERGDFNHERGERKVCRFKLEWLGNC GLNDETYGYKEGKPCIIKLNRLVLFKPKPPKNESLETYPVMKYNPVLPVQCT GKRDEDKDKIGNVEYFGLGNSPGFPLQYYPPYGGKLLQPKYLQPLLAVQFTNLT MDTEIRIECKAYGENIGYSEKDRFQGRFDVKIEVKS
Source	in vitro E.coli expression system
Target Names	ATP1B1
Protein Names	Recommended name: Sodium/potassium-transporting ATPase subunit beta-1 Alternative name(s): Sodium/potassium-dependent ATPase subunit beta-1
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.