



Recombinant *Torpedo californica* Sodium/potassium-transporting ATPase subunit beta-1 (atp1b1)

Product Code	CSB-CF002326TOT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P05029
Product Type	Transmembrane Protein
Immunogen Species	<i>Tetronarce californica</i> (Pacific electric ray) (<i>Torpedo californica</i>)
Sequence	MAREKSTDDGGGWKKFLWDSEKKQVLGRTGTSWFKIFVFLIFYGCLAGIFIG TIQVMLLTISDFEPKYQDRVAPPGLSHSPYAVKTEISFSVSNPNSYENHVNGLK ELLKNYNESKQDGNTPFEDCGVIPADYITRGPPIEESQGQKRVCERLLQWLKNC SGIDDPGYSEGGKPCIIAKLNRILGFYKPPKNGTDLPEALQANYNQYVLPPIHC QAKKEEDKVRIGTIEYFGMGGVGGFPLQYYPYKRLQKNYLQPLVGIQFTNL THNVELRVECKVFGDNIAYSEKDRSLGRFEVKIEVKS
Source	in vitro E.coli expression system
Target Names	atp1b1
Protein Names	Recommended name: Sodium/potassium-transporting ATPase subunit beta-1 Alternative name(s): Sodium/potassium-dependent ATPase subunit beta-1
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.