



Recombinant Pig Membrane cofactor protein (CD46)

Product Code	CSB-CF004939PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O02839
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	CDEPPKFESMRPQFLNTTYRPGDRVEYECRPGFQPMVPALPTFSVCQDDNT WSPLQEACRRKACSNLPDPLNGQVSYPNGDMLFGSKAQFTCNTGFYIIGAET VYCQVSGNVMWSEPSPLCEKILCKPPGEIPNGKYTNCHKDVFEYNEVVTYS CLSSTGPDEFSLVGEESLFCIGKDEWSSDPPECKVVKCPYPVVPNGEIVSGFG SKFYYKAEVVFKCNAGFTLHGRDTIVCGANSTWEPPEMPQCICKDSKPTDPPATP GPSHPGPPSPSDASPPKDAESLDGGIIAAIVGVLAIAVIAGGVYFFHHKYNKK RSK
Source	in vitro E.coli expression system
Target Names	CD46
Protein Names	Recommended name: Membrane cofactor protein Alternative name(s): CD_antigen= CD46
Expression Region	45-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.