



Recombinant Rat T-cell surface antigen CD2 (Cd2)

Product Code	CSB-CF004894RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P08921
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	RDSGTVWVGALGHGINLNIPNFQMTDDIDEVRWERGSTLVAEFKRKMKPFLKS GAFEILANGDLKIKNLTRDDSGTYNVTVYSTNGTRILDKALDLRILEMVSFKPMIY WECSNATLTCEVLEGTDVELKLYQGKEHLRSLRQKTMSYQWTNLRAPFKCKA VNRVSQESEMENVNCPKGLPLYLIVGVSAGLLLLVFFGALFIFCICKRKRNR RRKGEELEIKASRMSTVERGPKPHSTQASAPASQNPVASQAPPPGHHLQTP GHRPLPPSHRNREHQPKRPPPSGTQVHQKGPPLPRPRVQPKPPCGSGDV SLPPPN
Source	in vitro E.coli expression system
Target Names	Cd2
Protein Names	Recommended name: T-cell surface antigen CD2 Alternative name(s): LFA-2 LFA-3 receptor OX-34 antigen T-cell surface antigen T11/Leu-5 CD_antigen= CD2
Expression Region	23-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.