



Recombinant Saguinus oedipus Membrane cofactor protein (CD46)

Product Code	CSB-CF004939SWB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O62837
Product Type	Transmembrane Protein
Immunogen Species	Saguinus oedipus (Cotton-top tamarin)
Sequence	DACGPPPTFEAMELTSRPKPYKVGERVEYDCKKGYHHFAPFLTHSICDRNH TWLPISDEPCVRKVCHYIPNPLHGEAILANGSYSFGNQLHFICNDGYLIGKEIL YCELKGSDAVWVSGRPPICQKILCKPPPKNNGKHTFSDVDVFEYLDVAVTYS CDPAPGPDFSLIGESTIYCRDNSVWVSGDAPECKVVKCRFPVIENGKQIAGFGKK FYYKATVIFECDEGFHIGSDTVCSNSTWDPPVPCVKVSTSPATVSPTSSV PGYPNPDEGMLNSLDEWAIALIVIALVGVVAISFGLHRYLQRRKKKGGKADGTAE YATYQSKSATLAEQRS
Source	in vitro E.coli expression system
Target Names	CD46
Protein Names	Recommended name: Membrane cofactor protein Alternative name(s): CD_antigen= CD46
Expression Region	33-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.