



Recombinant Guinea pig Membrane cofactor protein (CD46)

Product Code	CSB-CF004939GU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70105
Product Type	Transmembrane Protein
Immunogen Species	Cavia porcellus (Guinea pig)
Sequence	CVLPPPFAMEPINPKPYEIGEKEVEYRCKKGYLRQPFYLMVATCEKNHSWVP ITDDGCIKKQCTYLNPPPKGRVEYINGTRTWGDIVHFSCVEGFYVSGIAALSCE LRGDNVDWNGRVPTCEKVLCSPPPKIQNGKYTFSDVQVFEYFEAVTYSCDAV QGPDKLSLVGNEVLYCAGHQKWSSAAPECKVVKCPLPVVKNQKQISGLGQTF FYQATVTFQCLPGFYFNGSSTVVCGSDNTWKPSIPECLKGPKPTHPTKPPVYN YPGYPNPREGIFDQELNLWIIILLIAVVGLALILLACRFFERKKKSPEKVMPAI QISDHMVNAPILMPSTVQNPMKR
Source	in vitro E.coli expression system
Target Names	CD46
Protein Names	Recommended name: Membrane cofactor protein Alternative name(s): CD_antigen= CD46
Expression Region	36-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.