



Recombinant Mouse Cytokine receptor common subunit gamma (Il2rg)

Product Code	CSB-CF011651MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34902
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	WSSKVLMSANEDIKADLILTSTAPEHLSAPTLPPEVQCFVFNIEYMNCTWNS SSEPQATNLTLHYRYKVSDNNTFQECSHYLFKSKEITSGCQIQKEDIQLYQTFVV QLQDPQKPQRRAVQKLNQNLVIPRAPENLTLSNLSESQLELRWKSRIKERC LQYLQYRSNRDRSWTELIVNHEPRFSLPSVDELKRYTFRVRSRYNPICGSSQ QWSKWSQPVHWGSHTVEENPSLFALEAVLIPVGTMGLIITLIFVYCWLERMPPPI PPIKNLEDLVTEYQGNFSAWSGVSKGLTESLQPDYSERFCHVSEIPPKGGALG EGPGGSPCSLHSPYWPPPCYSLKPEA
Source	in vitro E.coli expression system
Target Names	Il2rg
Protein Names	Recommended name: Cytokine receptor common subunit gamma Alternative name(s): Interleukin-2 receptor subunit gamma Short name= IL-2 receptor subunit gamma Short name= IL-2R subunit gamma Short name= IL-2RG gammaC p64 CD_antigen= CD132
Expression Region	23-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.