



# Recombinant Dog Cytokine receptor common subunit gamma (IL2RG)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-CF011651DO  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P40321  |
| <b>Product Type</b>      | Transmembrane Protein   |
| <b>Immunogen Species</b> | Canis lupus familiaris (Dog) (Canis familiaris)   |
| <b>Sequence</b>          | LNSTVPMPNGNEDITPDFFLTATPSETLSVSSLPLPEVQCFVFNVEYMNCTWN<br>SSSEPRPTNLTLHYWYKNSNDDKVQECGHYLSREVTAGCWLQKEEIHLYET<br>FVVQLRDPREPRRQSTQKLLQNLVIPWAPENLTLHNLSESQLELSWSNRHLD<br>HCLEHVQYRSDWDRSWTEQSVDRNSFSLPSVDGQKFYTFRVRSRYNPLC<br>GSAQRWSEWSHPIHWGSNTSKENPLFASEAVLIPLGSMGLIISLICVYYWLEERS<br>IPRIPTLKNLEDLVTEYHGNFSAWSGVSKGLAESLQPDYSEWLCHVSEIPPKG<br>GAPGEGPGGSPCSQHSPLYWAPPCYTLKPETGALIP |
| <b>Source</b>            | in vitro E.coli expression system   |
| <b>Target Names</b>      | IL2RG   |
| <b>Protein Names</b>     | Recommended name: Cytokine receptor common subunit gamma Alternative name(s): Interleukin-2 receptor subunit gamma Short name= IL-2 receptor subunit gamma Short name= IL-2R subunit gamma Short name= IL-2RG gammaC p64 CD_antigen= CD132  |
| <b>Expression Region</b> | 23-373  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |