



Recombinant Human Starch-binding domain-containing protein 1 (STBD1)

Product Code	CSB-CF022820HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95210
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MGAVWSALLVGGGLAGALFVWLLRGGPGDTGKDGDAEQEKDAPLGGAAIPG GHQSGSSGLSPGPSGQELVTKPEHLQESNGHLISKTKDLGKLQAASWRLQNP SREVCDNSREHVPSGQFPDTEAPATSETSNRSYSEVSRNESLESPMGEWG FQKGQEISAKAATCFAEKLPSNLLKNRAKEEMSLSDLNSQDRVDHEEWEMV PRHSSWGDVGVGGSLKAPVLNLNQGMDNGRSTLVEARGQQVHGKMERVAV MPAGSQQVSVRFQVHYVTSTDVQFIAVTGDHECLGRWNTYIPLHYNKDGFW HSIFLPADTVVEWKFVLENGGVTRWEECSNRFLETGHEDKVVHAWWGIH
Source	in vitro E.coli expression system
Target Names	STBD1
Protein Names	Recommended name: Starch-binding domain-containing protein 1 Alternative name(s): Genethonin-1
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.