



Recombinant Mouse Neuroplastin (Nptn)

Product Code	CSB-CF016028MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P97300
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	QNAGFVKSPMSETKLTGDAFELYCDVVGSPTEIQWWYAEVNRAESFRQLW DGARKRRVTNTAYGSNGVSVLRITRLTLEDSTYECRASNDPKRNDLRQNP SITWIRAQATISVLQKPRIVTSEEVIRESLLPVTLCNLTSSHTLMYSYWTRNG VELTATRKNASNMEYRINKPRAEDSGEYHCVYHFVSAPKANATIEVKAAPDITG HKRSENKNEGQDAMMYCKSVGYPHPEWIWRKKENGVFEEISNSSGRFFITNK ENYTELSIVNLQITEDPGEYECNATNSIGSASVSTVLRVRSHLAPLWPFLGILAEI IILVVIIVVYEKRKRKPDEVPDDDEPAGPMKTNSTNNHKDKNLRQRNTN
Source	in vitro E.coli expression system
Target Names	Nptn
Protein Names	Recommended name: Neuroplastin Alternative name(s): Stromal cell-derived receptor 1 Short name= SDR-1
Expression Region	29-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.