



Recombinant Human Uncharacterized protein C21orf63 (C21orf63)

Product Code	CSB-CF003777HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P58658
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	LTDFSGYLTKLLQNHTTYACDGYLNLQCPRHSTISVQSAFYGQDYQMCSSQ KPASQREDSLTCVAATTFQKVLDECQNQRACHLLVNSRVFGPDLCPGSSKYL LVSFKCQPNELKNKTVCEDQELKLHCHESKFLNIYSATYGRRTQERDICSSKA ERLPPFDCLSYSALQVLSRRCYGKQRCKIIVNNHHFGSPCLPGVKKYLTVTYA CVPKNILTAIDPAIANLKP SLKQKDGEYGINFDPSGSKVLRKDGILVSNSLAFA YIRAHPERAALLFVSSVCIGLALTLCALVIRESCAKDFRDLQLGREQLVPGSDKV EEDSEDEEEEEEDPSESDFPGELSGFCRTSYPIYSSIEAAELAERIERREQUIQEI WMNSGLDTS LPRNMGQFY
Source	in vitro E.coli expression system
Target Names	EVA1C
Protein Names	Recommended name: Uncharacterized protein C21orf63 Alternative name(s): SUE21
Expression Region	49-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.