



Recombinant Mouse Interleukin-5 receptor subunit alpha (Il5ra)

Product Code	CSB-CF011663MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P21183
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	DLLNHKKFLLLPPVNFTIKATGLAQVLLHWDPNPDQEQRHVDLEYHVKINAPQE DEYDTRKTESKCVTPLHEGFAASVRTILKSSHTTLASSWWSAELKAPPGSPGT SVTNLTCTTHTVVSSHTHLRPYQVSLRCTLVVGKDAPEDTQYFLYYRFGVLTE KCQEYSRDALNRNTACWFPRTFINSKGFQQLAVHINGSSKRAAIKPFDQLFSP L AIDQVNPPRNVTVEIESNSLYIQWEKPLSAFPDHCFFNYELKIYNTKNGHIQKEKL IANKFISKIDDVSTYSIQVRAAVSSPCRMPPGRWGEWSQPIYVGKERKSLVEWH LIVLPTAACFVLLIFSLICRVCHLWTRLFPPVPAPKSNIKDLPVVTEYEKPSNETKI EVVHCVEEVGFVVMGNSTF
Source	in vitro E.coli expression system
Target Names	Il5ra
Protein Names	Recommended name: Interleukin-5 receptor subunit alpha Short name= IL-5 receptor subunit alpha Short name= IL-5R subunit alpha Short name= IL-5R- alpha Short name= IL-5RA Alternative name(s): CD_antigen= CD125
Expression Region	18-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.