



Recombinant Mouse Interleukin-13 receptor subunit alpha-1 (Il13ra1)

Product Code	CSB-CF011591MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O09030
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	ATEVQPPVTNLSVSVENLCTIITWSPPEGASPNCTLRYSHFDDQQDKKIAP ETHRKEELPLDEKICLQVGSQCSANESEKPSPLVKKCISPPEGDPESAVTELKC IWHNLSYMKCSWLPGRNTSPDTHYTLYYWYSSLEKSRQCENIYREGQHIACS FKLTKVEPSFEHQNVQIMVKDNAGKIRPCKIVSLTSYVKPDPPHIKHLKNGA LLVQWKNPQNFRSRCLTYEVEVNNTQTDRHNILEVEEDKCQNSESDRNMEGT SCFQLPGVLADAVYTVRVRVKTNLKCFDDNKLWSDWSEAQSIGKEQNSTFYT TMLLTIPVFVAVAVIILLFYLKRLKIIFFPIPDGPKIFKEMFGDQNDLHWKKYDI YEKQSKEETDSVVLIENTLKKAAP
Source	in vitro E.coli expression system
Target Names	Il13ra1
Protein Names	Recommended name: Interleukin-13 receptor subunit alpha-1 Short name= IL-13 receptor subunit alpha-1 Short name= IL-13R subunit alpha-1 Short name= IL-13R-alpha-1 Short name= IL-13RA1 Alternative name(s): Interleukin-13-binding protein N
Expression Region	26-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.