



Recombinant Chicken Tapasin (TAPBP)

Product Code	CSB-CF023122CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O73895
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	QFRVEDAASPPPPPAPVRCALLEGVGRGGGLPGGGNARPALLRFGGDAETP PEPGPEPEVTFNVSDPWGTLTPLGVPPRTPPSCELNPTNPQTGSDPWSRPLH PDARSPPTAGGQWWVAAVGTPQYGV TALLQGGMGTEGTITAAVALAVLTHTP TLRARVGSPIHLHCAFAAPPSSFVLEWRHQNRGAGRVLLAYDSSTARAPRAH PGAELLGLTRDGDGVTAVTLRLARPSPGDEGTYICSVFLPHGHTQTVLQLHVF EPPKVTLSPKNLVVAPGTS AELRCHVSGFYPLDVTVTWQRRAGGSGTSQSPR DTVMDSWTSGHRQAADGTYSRTAAARLIPARPQHGGDIYSCVVHTALAKPM RVSVRLLLAGTEGPHLEDITGLFLVAFVLCGLIRWLYPKAAPKEETKKSQ
Source	in vitro E.coli expression system
Target Names	TAPBP
Protein Names	Recommended name: Tapasin Short name= TPN Short name= TPSN Alternative name(s): TAP-associated protein TAP-binding protein
Expression Region	16-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.