



Recombinant Danio rerio Heparan-sulfate 6-O-sulfotransferase 3-B (hs6st3b)

Product Code	CSB-CF010759DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0MGZ7
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Transmembrane Proteins
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MNDKPNKWIFIPILAILFVMIGYQYVCPAGGQACHFRTGDKLVRIAPLATPDPTT DDLYREQDPEEDNPPKCASKFNFTERDLTRDVFDFNIKDDVIVFLHIQKTGGTT FGRHLVRNIRLEQPCDCKAGQKKCTCHRPQKQESWLFSRFSTGWSCGLHAD WTELTNCVPVIMDKRQPPKRKRNFYYITMLRDPVSRYLSEWKHVQRGATWKT SLHMCDGRSPTQDELPTCYNGDDWSGVTLHDFMDCPSNLANNRQVRMLADL SLVGCYNLSTMNESERNPILLASAKSNLKNMAFYGLTEFQRKTQYLFERTFHL RFISAFTQINSTRAANVELRDDMRSRIEQLNMLDMQLYEFADLFLQRYQFVR QRERQEERLKRREERRWIRERRVNQSKEPIVENQTRVTTTETYASQVVRW
Source	in vitro E.coli expression system
Target Names	hs6st3b
Protein Names	Recommended name: Heparan-sulfate 6-O-sulfotransferase 3 Short name= HS 6-OST-3 EC= 2.8.2.-
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.