



# Recombinant Human Killer cell immunoglobulin-like receptor 3DL2 (KIR3DL2)

<b>Product Code</b>	CSB-CF012365HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P43630
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	LMGGQDKPFLSARPSTVVPRGGHVALQCHYRRGFNNFMLYKEDRSHVPIFHG RIFQESFIMGPVTPAHAGTYRCRGSRPHSLTGWSAPSNPLVIMVTGNHRKPSL LAHPGPLLKSGETVILQCWSDVMFEHFFLHREGISEDPSRLVGQIHDGVSKAN FSIGPLMPVLAGTYRCYGSVPHSPYQLSAPSDPLDIVITGLYEKPSLSAQPGPT VQAGENVTLSCSSWSSYDIYHLSREGEAHERRLRVAVPKVNRTFQADFPLGPA THGGTYRCFGSFRALPCVWSNSSDPLLVSVTGNPSSSWPSPTPESSKSGICR HLHVLIGTSVVIFLIFILLFFLLYRWCSNKKNAAVMDQEPAGDRTVNRQDSDEQ DPQEVTYAQLDHCVFIQRKISRPSQRPKTPLTDTSVYTELPNAEPRSKVVSCP RAPQSGLEGVF
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	KIR3DL2
<b>Protein Names</b>	Recommended name: Killer cell immunoglobulin-like receptor 3DL2 Alternative name(s): CD158 antigen-like family member K MHC class I NK cell receptor Natural killer-associated transcript 4 Short name= NKAT-4 p70 natural killer cell receptor cl
<b>Expression Region</b>	22-455
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.