



Recombinant Synaptotagmin-1 (snt-1)

Product Code	CSB-CF023029CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34693
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MVKLDFSSQDEENDEDLTKEFVRDEAPMEETTSEAVKQIATTTKETLKDVVVN KVIDVKDVVKEKVMQQTGMPEWAFVFLGFVFILLVLACAFCLIRKLFGKKRHGE KNKKGGLKGFFGKGQDVVDGKNIQGMAQDLEELGDAMEQNEKEQAEKEEV KLGRIQYKLDYDFQQGQLTVTVIQAEDLPGMDMSGTSDPYVKLYLLPEKKKKV ETKVHRKTLNPVFNETFIFKVAFNEITAKTLVFAIYDFDRFSKHDQIGQVLIPLGKI DLGAVIEEWKDIAPPPDDKEAEKSLGDICFSLRYVPTAGKLTVVILEAKNLKKMD VGGLSDPYVKIVLMQGGKRLKKKTSIKKCTLNPYYNESFSFEVPFEQIQKVSL MITVMDYDKLGSNDAIGRCLLGCGNGTGAELRHWMDMLASPRRPIAQWHTLGP VEEEGDKKDDKK
Source	in vitro E.coli expression system
Target Names	snt-1
Protein Names	Recommended name: Synaptotagmin-1 Alternative name(s): Synaptotagmin I
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.