



Recombinant Human Syndecan-3 (SDC3)

Product Code	CSB-CF020890HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75056
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MKPGPPHRAGAAHGAGAGAGAAAGPGARGLLLPPLLLLLLAGRAAGAQRWR SENERPVDLEGSGDDDSFPDDELDDLYSGSGSGYFEQESGIETAMRFSPDV ALAVSTTPAVLPTTNIQPVGTPFEELPSERPTLEPATSPLVVTEVPEEPSQRAT TVSTTMATTAATSTGDPTVATVPATVATATPSTPAAPPFTATTAVIRTTGVRRL LPLPLTTVATARATTPEAPSPPTTAAVLDTAAPTPLRVSTATSRPRALPRPATT QEPDIPERSTLPLGTTAPGPTVAQTPTPETFLTTRDEPEVPVSGGPGSDFEL PEEETTQPDTANEVVAVGGAAAKASSPPGTLPKGARPGPGLLDNAIDSGSSA AQLPQKSILERKEVLVAVIVGGVVGALFAAFLVTLIIYRMKKKDEGSYTLLEPKQ ASVTYQKPKDKQEEFYA
Source	in vitro E.coli expression system
Target Names	SDC3
Protein Names	Recommended name: Syndecan-3 Short name= SYND3
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.