



Recombinant Pig Interleukin-6 receptor subunit alpha (IL6R)

Product Code	CSB-CF011665PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O18796
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	LAPRGCSKLEVAQDVLTSPLGASVTLTCTPGGEPGDNATIHWVLRNQVTGSPD GRPAGVGRLLLLKSVQLSDSGNYSCYQDGVPAAGSVRLLVDAPPEEPQLSCFR KSPLSNVGCEWRPRSPPTTKAVLLVRKFQNSPVEDFQEPQCQYSLEAQRFF CQLAVPEGDNSFHIVTLCVANSAGSQSSTPQTFEGYGILQPDPPVNITVSAVD RNPRWLSVTWQDPPSWNSYFYRLQFELRYRAERSKTFTTWMVKELQHHCIIH DAWSGMRHVVQLRAQEEFGHGLWSEWSQEVGTGIPWTESRSSPAETELPLST QAPTTNEDDEDISSKESANATSLPVQDSASVPLPTFLVAGGSLAFGTLLCIGIIL RFKKTGQLQALKEGKTNMHPYPYSLGQLVPERPKSTPVLVPLISPPVSPNSLGD NTSRNSRPEARQPQSPYDVSNRDYFFPR
Source	in vitro E.coli expression system
Target Names	IL6R
Protein Names	Recommended name: Interleukin-6 receptor subunit alpha Short name= IL-6 receptor subunit alpha Short name= IL-6R subunit alpha Short name= IL-6R-alpha Short name= IL-6RA Alternative name(s): IL-6R 1 CD_antigen= CD126
Expression Region	20-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.