



# Recombinant Chlorocebus aethiops Hepatitis A virus cellular receptor 1 (HAVCR1)

<b>Product Code</b>	CSB-CF010144DSU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O46598
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
<b>Sequence</b>	<p>DSVNVDGVAGLSITLPCRYNGAITS MCWNRGTCSVFSCP DGIVWTNGTHV TY            RKETRYKLLGNLSRRDVS LTIANTAVSDSGIYCCR VKHSGWFNDMKITIS LKIGP            PRVTIPIVRTVRTSTTVPTTTTTLPTTTTLPTTTTLPTTMTLPTTTTLPMTTTLPT            TTTVPMTTTLPTTLPTTTTLPTTLPTTTTLPTTLPTTTTLPTTMTLPMTTTLPTTT            TLPTTTTLPTTTTLPTTTTLPTTLPTTMTLPTTTTLPTTMTLPMTTTLPTTTTLPTT            TTLPTTTMVFSTFVPPTPLPMQDHEPVATSPSSAQPAETHPV TLLGATRTRQPTS            SPLYSYTTDGS DTVTESSDGLWNNNQ TQLSPEHSPQMVNTTEGIYAGVCISVL            VLLAVLGVVIAKKYFFKKEIQQLSVSFSNHQFKTLQNAV KKEVHAEDNIYIENNL            YAMNQDPV VLFESLRP</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	HAVCR1
<b>Protein Names</b>	Recommended name: Hepatitis A virus cellular receptor 1 Short name= HAVcr-1 Alternative name(s): T-cell immunoglobulin and mucin domain-containing protein 1 Short name= TIMD-1 T-cell membrane protein 1 TIM-1
<b>Expression Region</b>	25-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.