



Recombinant Mouse Interferon gamma receptor 1 (Ifngr1)

Product Code	CSB-CF011051MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15261
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	<p>ALTSTEDPEPPSVPVPTNVLIKSYNLNPVVCWEYQNMSQTPIFTVQVKVYSGS WTDSCNTN ISDHCCNIYEQIMYPDVSARVVKAKVGQKESDYARSKEFLMCLKGKVGPPG LEIRRKKE EQLSVLVFHPVVVNGESQGTMTFGDGSTCYTFDYTVYVEHNRSGEILHTKHT VEKEECNE TLCELNISVSTLDSRYCISVDGISSFQVVRTEKSKDVCIPPFHDDRKDSIWILVV APLTV FTVVILVFAYWYTKKNSFKRKSIMLPKSLLSVVKSATLETKPESKYSLVTPHQPA VLESE TVICEEPLSTVTAPDSPEAAEQEELSKETKALEAGGSTSAMTPDSPPTPTQRR SFSLLSS NQSGPCSLTAYHSRNGSDSGLVVGSGSSISDLESLPNNNSETKMAEHDPPPVR KAPMASGY DKPHMLVDVLVDVGGKESLMGYRLTGAEQELS</p>
Source	in vitro E.coli expression system
Target Names	Ifngr1
Protein Names	Recommended name: Interferon gamma receptor 1 Short name= IFN-gamma receptor 1 Short name= IFN-gamma-R1 Alternative name(s): CD_antigen= CD119
Expression Region	26-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.