



Recombinant Human Interferon alpha/beta receptor 2 (IFNAR2)

Product Code	CSB-CF011047HU(A4)
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48551
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	ISYDSPDYTDESCTFKISLRNFRSILSWELKNHSIVPTHYTLTYTIMSKPEDLKVV KNCANTTRSFCDLTDEWRSTHEAYVTVLEGFSGNTTLFSCSHNFWLAIDMSFE PPEFEIVGFTNHINVMVKFPSIVEEELQFDLSLVIEEQSEGIVKKHKPEIKGNMS GNFTYIIDKLIPNTNYCVSVYLEHSDEQAVIKSPLKCTLLPPGQESAEAKIG GIITVFLIALVLTSTIVTLKWIGYICLRNSLPKVLNFHNFLAWPFPNLPPLEAMD VEVIIYNRKKKWDYNYDDESDDTEAAPRTSGGGYTMHGLTVRPLGQASAT STESQLIDPESEEEPDLPEDVDLPTMPKDSPQQLELLSGPCERRKSPLQDPF PEEDYSSTEGSGGRITFNVDLNSVFLRVLDDESDDLEAPLMLSSHLEEMVDP EDPDNVQSNHLLASGEGTQPTFPSPSSEGLWSEDAPSDQSDTSESVDLGD GYIMR
Source	in vitro E.coli expression system
Target Names	IFNAR2
Protein Names	Recommended name: Interferon alpha/beta receptor 2 Short name= IFN-R-2 Short name= IFN-alpha binding protein Short name= IFN-alpha/beta receptor 2Alternative name(s): Interferon alpha binding protein Type I interferon receptor 2
Expression Region	27-515
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.