



# Recombinant Mouse Activin receptor type-2A (Acvr2a)

<b>Product Code</b>	CSB-CF001260MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P27038
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Sequence</b>	<p>AILGRSETQECLFFNANWERDRTNQTGVEPCYGDKDKRRHCFATWKNISGSI  EIVKQGCWLDDINCYDRDTCIEKKDSPEVYFCCCEGNCNEKFSYFPEMEVT  QPTSNPVTPKPPYYNILLSLVPLMLIAGIVICAFWVYRHHKMAYPPVLVPTQDP  GPPPPSPLLGLKPLQLLEVKARGRFVCVWKAQLLNEYVAVKIFPIQDKQSWQN  EYEVYSLPGMKHENILQFIGAEKRGTSVDVDLWLITAFHEKGSLSDFLKANVVS  WNECHIAETMARGLAYLHEDIPGLKDGHKPAISHRDIKSKNVLLKNNLTACIAD  FGLALKFEAGKSAGDTHGQVGTRRYMAPEVLEGAINFQRDAFLRIDMYAMGL  VLWELASRCTAADGPVDEYMLPFEEEEIGQHPSLEDMQEVVVHKKKRPVLRDY  WQKHAGMAMLCETIEECWDHDAEARLSAGCVGERITQMQRRLTNIITTEDIVTV  VTMVTNVDFPPKESSL</p>
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Acvr2a
<b>Protein Names</b>	Recommended name: Activin receptor type-2A EC= 2.7.11.30 Alternative name(s): Activin receptor type IIA Short name= ACTR-IIA
<b>Expression Region</b>	20-513
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.