



Recombinant Human Syntaphilin (SNPH)

Product Code	CSB-CF022317HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15079
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAMSLPGSRRTSAGSRRRTSPPVSVRDAYGTSSLSSSSNSGSGSYKGSDDSSPT PRRSMKYTLCSDNHGIKPPTPEQYLTPQQKEVCIRHLKARLKDTQDRLQDRD TEIDDLKTQLSRMQEDWIEEECHRVEAQLALKEARKEIKQLKQVIDTVKNNLIDK DKGLQKYFVDINIQNKKLETLLHSMEVAQNGMAKEDGTGESAGGSPARSLTR SSTYTKLSDPAVCGDRQPGDPSSGSAEDGADSGFAAADDTLSRTDALEASSL LSSGVDCGTEETSLHSSFGLGPRFPASNTYEKLLCGMEAGVQASCMQERAIQ TDFVQYQPDLDTILEKVTQAQVCGTDPESGDRCPELDAHPSGPRDPNSAVVV TVGDELEAPEPITRGPTPQRPGANPNPGQSVSVVCPMEEEEEAAVAEKEPKS YWSRHYIVDLLAVVPAVPTVAWLRCRSQRRRQQQPIYNISLLRGCCCTVALHSIR RISCRSLSQPSPSPAGGGSQ
Source	in vitro E.coli expression system
Target Names	SNPH
Protein Names	Recommended name: Syntaphilin
Expression Region	1-494
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.