



Recombinant Rat Interleukin-2 receptor subunit beta (Il2rb)

Product Code	CSB-CF011650RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P26896
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	AVNDCSHLKCFYNSRANVSCMWSPEEALNVTSCHIHAKSDMRHWNKTCELT PVRQASWACNLILGPLPDSQSLTSVDLLSLSVVCWEEKGWRRVKTCTFHPFD NLRLIAPHSLQVLHIETRRCNISWEVSQVSHYVNPYLEFEARRLLDRSWEDA SVFSLKQRQQWIFLETLPDTSYELQVRVIAQRGKTRTWSPWSQPMAFRTRP ADPKEIFPLPWLRCLLLVLGCFFGFLSCVCVLVKCRYLGPWLKTLKCHIPDPS EFFSQLSSQHGGDLQKWLSSPVPQSFPSPTGSAPEISPLEVDRDSDTKMQL LFQKEKASSPSPSGHSQASCFTNQGYFFFHLSNALEIESCQVYFTYDPCMEED VEEDGPRLPEESPLPPLLPTGEQDDYCAFPPRDDLLLFSPSMSTPNTAYGNS ITPEERPPLSLQEGLPSLSPDLMLGLQHLELELGDDGEGMSTNSSGQQASV PEAALMGTTKTEARPVLTNTDAYLSLQELQAQDSAHLI
Source	in vitro E.coli expression system
Target Names	Il2rb
Protein Names	Recommended name: Interleukin-2 receptor subunit beta Short name= IL-2 receptor subunit beta Short name= IL-2R subunit beta Short name= IL-2RB Alternative name(s): High affinity IL-2 receptor subunit beta p70-75 CD_antigen= CD122
Expression Region	27-537
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.