



Recombinant Pelophylax nigromaculatus Tyrosinase (TYR)

Product Code	CSB-CF025394PEM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q04604
Product Type	Transmembrane Protein
Immunogen Species	Pelophylax nigromaculatus (Black-spotted frog) (Rana nigromaculata)
Sequence	GIDGHRGGRGTHQSNIRIYSDSAPPWFTKEDISAMRFLSDSRIGHIKQNLLLFE SDQTPL
Source	in vitro E.coli expression system
Target Names	TYR
Protein Names	Recommended name: Tyrosinase EC= 1.14.18.1 Alternative name(s): Monophenol monooxygenase
Expression Region	20-532
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.