



Recombinant Mouse Interleukin-1 receptor-like 1 (Il1rl1)

Product Code	CSB-CF011626MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P14719
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	SKSSWGLENEALIVRCPQRGRSTYPVEWYYSDTNESIPTQKRNRIFVSRDLK FLPARVEDSGIYACVIRSPNLNKTGYLNVTIHKKPPSCNIPDYLMYSTVRGSDK NFKITCPTIDLYNWTAPVQWFKNCKALQEPRFRAHRSYLFIDNVTHDDEGDYT CQFTHAENGNTNYIVTATRSFTVEEKGFMSFPVITNPPYNHTMEVEIGKPASIA SACFGKGSFLADVLWQINKTVVGNFGEARIQEEGRNESSNDMDCLTSVL RITGVTEKDLSLEYDCLALNLHGMIRHTIRLRKQPIDHRSIYYIVAGCSLLMFI NVLVIVLKVFWIEVALFWRDIVTPYKTRNDGKLYDAYIIPRVFRGSAAGTHSVE YFVHHTLPDVLENKCGYKLCIYGRDLLPGQDAATVVESSIQNSRRQVFLAPH MMHSKEFAYEQEIALHSALIQNNSKVIEMEPLEASRLQVGDQLQDSLQHLVKI QGTIKWREDHVADKQSLSSKFWKHVRYQMPVPERASKTASVAAPLSGKACLD LKHF
Source	in vitro E.coli expression system
Target Names	Il1rl1
Protein Names	Recommended name: Interleukin-1 receptor-like 1 Alternative name(s): Lymphocyte antigen 84 Protein ST2 Protein T1
Expression Region	27-567
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.