



Recombinant Macaca mulatta (Rhesus macaque) Interleukin-1 receptor accessory protein (IL1RAP)

Product Code	CSB-CF011623MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59822
Product Type	Transmembrane Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Sequence	SERCDDWGLDTMRQIQVFEDEPARIKCPLFEHFLKFNYSTAHSAGLTLIWYWT RQDRDLEEPINFRLPENRISKEKDVLWFRPTLLNDTGNYTCMLRNTTYCSKVA FPLEVVQKDSCFNSPMKLPVHKLYIEYGIQRITCPNVDGYFPSSVKPTITWYMG CYKIQNFNNVIPEGMNL SFLIAFISNNGNYTCVVTPENGRTFHLTRTLTVKVVG SPKNAVPPVIHSPNDHVVEKEPGEELLIPCTVYFSFLMDSRNEVWWTIDGKK PDDIPIDVTINESISHSRTEDETRTQILSIKKVTSEDLKRSYVCHARSAKGEVAKA ATVKQKVPAPRYTVELACGFGATVLLVVILIVVYHVYWLEMVLFYRAHFGTDET ILDGKEYDIYVSYARNAEEEEFVLLTLRGVLENEFGYKLCIFDRDSLPGGIVTDE TLSFIQKSRLLVVLSPNYVLQGTQALLELKAGLENMASQGNINVILVQYKAVK ETKVKELKRAKTVLTVIKWKGEKSKYPQGRFWKQLQVAMPVKKSPRRSSSDE QGLSYSSLKNV
Source	in vitro E.coli expression system
Target Names	IL1RAP
Protein Names	Recommended name: Interleukin-1 receptor accessory protein Short name= IL-1 receptor accessory protein Short name= IL-1RAcP
Expression Region	21-570
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.