



Recombinant Chicken CD166 antigen (ALCAM)

Product Code	CSB-CF001562CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42292
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	LYTVNAVYGDITIMPCRLEVPDGLMFGKWKYEMPNGSPVFIAFRSSTKKNVQ YDDVPDYKDRLSLSENYTLSEIKNARISDEKRFVCMLVTEDDVSEEPTVVKVFKQ PSQPEILHQADFLETEKLKMLGECVVRDSYPEGNVTWYKNGRVLQPVEEVVVI NLRKVENRSTGLFTMTSSLQYMPKEDANAKFTCIVTYHGPSGQKTIQSEPVV FDVHYPTEKVTIRVLSQSSTIKEGDNVTLKCSGNGNPPPPQEFLFYIPGETEGIR SSDTYVMTDVRRNATGEYKCSLIDKSMMDATTITVHYLDLQLTPSGEVTKQIGE ALPVSCTISSRNATVFWIKDNTRMKTSPSFSSLQYQDAGNYICETTLQEVEGL KKRKTLLIVEGKPKQIKMTKKTNTNKMSKTIVCHVEGFPAVQWTVTGSGSLI NKTEETKYVNGKFSSKIIIAPEENVTLTICIAENELERTVTSLNVS AISIPEYDEPED RNDNSEKVNDAQKLIVGIVVGLLLVALVAGVVYWLIVKKSKTASKHVVDKDLG NIEENKLEENNHKSET
Source	in vitro E.coli expression system
Target Names	ALCAM
Protein Names	Recommended name: CD166 antigen Alternative name(s): Activated leukocyte cell adhesion molecule BEN glycoprotein Protein DM-GRASP Protein JC7 SC1 glycoprotein CD_antigen= CD166
Expression Region	34-588
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.