



Recombinant Mouse Mucin-13 (Muc13)

Product Code	CSB-CF015217MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P19467
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	SSSQASSTTSSSGGTSPPTTVQSQSPGSSSQASTTTSSSGGASPPTTVQSQS PGSSSQASTTTSSSGGASPPTTVQSQSPGSSSQASTTTSSSGGASPPTTVQSQ QSPGSSSQASTTTSSSGGASPPTTVQSQSPGSSSQASTTTSSSGGASPPTTV QSQSPGSSSQVSTTTSSSGGASPPTTVQSQSPGSSSQPGPTQPSGGASSST VPSGGSTGPS DLCNPNPCKGTASCVKLHSHKHFCLCLEGYYNSSLSSCVKGT TFPGDISMSVSETANLEDENSVGYQELYNSVTDFETFNKTQDYGQTVIIKVST APRSARSAMRDATKDVSVSVNIFGADTKETEKSVSSAIETAIKTSGNVKDYV SINLCDHYGCVGNDSSKQCQDILQCTCKPGLDRLNPQVPCVAVTCSQPCNAE EKEQCLKMDNGVMDCVCMPCGYQRANGNRKCEECFPGYSGMNCCKDQFQLIL TIVGTIAGALILILLIAFIVSARSKNKKKDGEEQRLIEDDFHNLRLRQTGFSNLGAD NSIFPKVRTGVPSQTPNPYANQRSMPRPDY
Source	in vitro E.coli expression system
Target Names	Muc13
Protein Names	Recommended name: Mucin-13 Short name= MUC-13 Alternative name(s): Cell surface antigen 114/A10 Lymphocyte antigen 64
Expression Region	18-573
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.