



Recombinant Mouse Delta-like protein 3 (Dll3)

Product Code	CSB-CF006948MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O88516
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	HSFGPGPGLGTPRSPCNARGPCRLFFRVCLKPGVSQEATESLCALGAALSTS VPVYTEHPGESAAALPLPDGLVRVPFRDAWPGTFSLVETWREQLGEHAGGP AWNLLARVVGRRRRLAAGGPWARDVQRTGTWELHFSYRARCEPPAVGAACA RLCRSRSAPSRCGPGLRPCTPFPECEAPSVCRPGCSPEHGYCEEPDECRC LEGWTGPLCTVPVSTSSCLNSRVPGPASTGCLLPGPGPCDGNPCANGGSCS ETSGSFECACPRGFYGLRCEVSGVTCADGPCFNGGLCVGGEDPDSAYVCHC PPGFQGSNCEKRVDRCSLQPCQNGGLCLDLGHALRCRCRAGFAGPRCEHDL DDCAGRACANGGTCVEGGGSRRCSALGFGGDRDCRERADPCASRPCAHGG RCYAHFSGLVACAPGYMGVRCFAVRPDGADAVPAAPRGLRQADPQRFL PPALGLLVAAGLAGAALLVIHVRRRGPGQDTGTRLLSGTREPSVHTLPDALNN LRLQDGAGDGPSSADWNHPEDGDSRSIYVIPAPSIYAREDWLIQVLF
Source	in vitro E.coli expression system
Target Names	Dll3
Protein Names	Recommended name: Delta-like protein 3 Alternative name(s): Drosophila Delta homolog 3 Short name= Delta3 Short name= M-Delta-3
Expression Region	33-592
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.