



Recombinant Human Interleukin-18 receptor accessory protein (IL18RAP)

Product Code	CSB-CF011611HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95256
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	<p>FNISGCSTKKLLWTYSTRSEEEFVLFCDLPEPQKSHFCHRNRLSPKQVPEHLP FMGSNDLSDVQWYQQPSNGDPLEDIRKSYPHIIQDKCTLHFLTPGVNNSGSYI CRPKMIKSPYDVACCVKMILEVKPQTNASCEYSASHKQDLLLGSTGSGISCP CQSDAQSPAVTWYKNGKLLSVERSNRIVVDEVYDYHQGTVCYDYSQSDTVSS WTVRAVVQVRTIVGDTKLKPDILDPEVDTLEVELGKPLTISCKARFGFERVFN PVIKWIYKDSLEWEVSVPEAKSIKSTLKDEIIERNIILEKVTQRDLRRKFCFVQ N SIGNTTQSVQLKEKRGVLLYILLGTIGTLVAVLAASALLYRHWIEIVLLYR TYQS KDQTLGDKKDFDAFVSYAKWSSFPSEATSSLSEEHLALSFPDVLENKY GYSL CLLERDVAPGGVYAEDIVSIIKRSRRGIFILSPNYVNGPSIFELQAAV NLALDDQT LKLILIKFCYFQEPESLPHLVKKALRVLPTVTWRGLKSVPPNS RFWAKMRYHM PVKNSQGFTWNQLRITSRIFQWKGLSRTETTGRSSQPKEW</p>
Source	in vitro E.coli expression system
Target Names	IL18RAP
Protein Names	Recommended name: Interleukin-18 receptor accessory protein Short name= IL-18 receptor accessory protein Short name= IL-18RAcP Alternative name(s): Accessory protein-like Short name= AcPL CD218 antigen-like family member B CDw218b IL-1R
Expression Region	20-599
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.