



Recombinant Human E-selectin (SELE)

| | |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-CF020975HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P16581 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Homo sapiens (Human) |
| Sequence | WSYNTSTEAMTYDEASAYCQQRYTHLVAIQNKEEIEYLNSILSYSPSYWIGIR KVNNVVWVWGTQKPLTEEAKNWAPGEPNNRQKDEDCVEIYIKREKDVGMWN DERCSKKLALCYTAACNTNTSCSGHGECVETINNYTCKCDPGFSGLKCEQIVN CTALESPEHGSLVCSHPLGNFSYNSSCSISCDRGLPSSMETMQCMSSGEW SAIPACNVVECDAVTNPANGFVECFQNPGSFPWNTTCTFDCEEFGFELMGAQ SLQCTSSGNWDNEKPTCKAVTCRAVRQPQNGSVRCSHSPAGEFTFKSSCNF TCEEFGFMLQGPAQVECTTQGQWTQQIPVCEAFQCTALSNPERGYMNCLPSA SGSFRYGSSCEFSCEQGFLKGSKRLQCGPTGEWDNEKPTCEAVRCDAVHQ PPKGLVRCACHSPIGEFTYKSSCAFSCEEFGFELHGSTQLECTSQGQWTEEVPS CQVVKCSSLAVPGKINMSCSGEPVFGTVCKFACPEGWTLNGSAARTCGATG HWSGLLPTCEAPTESNIPLVAGLSAAGLSLLTLAPFLLWLRKCLRKAKKFVPAS SCQSLESDGSYQKPSYIL |
| Source | in vitro E.coli expression system |
| Target Names | SELE |
| Protein Names | Recommended name: E-selectin Alternative name(s): CD62 antigen-like family member E Endothelial leukocyte adhesion molecule 1 Short name= ELAM-1 Leukocyte-endothelial cell adhesion molecule 2 Short name= LECAM2 CD_antigen= CD62E |
| Expression Region | 22-610 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |