



Recombinant Rabbit Prolactin receptor (PRLR)

Product Code	CSB-CF018727RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P14787
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	QSPPGKPFIFKCRSPEKETFTCWWRPGADGGLPTNYTLTYHKEGETITHECPD YKTGGPNSCYFSKKHTSIWTIYIITVNATNQMGSSVSDPRYVDVTYIVEPDPV NLTLEVKHPEDRKPYLWVKWLPPTLVDVRSWLTQYEIRLKPEKAAEWETHF AGQQTQFKILSLYPGQKYLQVRCKPDHGFWSVWSPESSIQIPNDFTMKDITV WIFVAVLSTIICLIMVWAVALKGYSMVTCIFPPVPGPKIKGFDTHLLEKKGKSEELL SAFGCQDFPPTADCEDLLVEFLEVDDSEDQQLMPAHSKEHSGPGMKPTDLDP DNDSGRGSCDSPSLSEKCEEPQANPSTFHTPEVIEQPEKPKANVTHTWDPQ TISLVGKMPYLSVNGSKSSTWPLLQPGQHNTNSPYHNIADMCKLATS�DKIDK DALQSSKTTEAAGEEKATKQREVESSHKAEQDTGWLLPKEKPPFISPKPLDY VEIHKNKDGALSLLLKQKENGDQTGKAGTPETSKEYAKVSRVMDNNILVLVQ DPGAQNVALFEESTKEAPSPSQNQAEDLSSFSTAPSDCRLQQGGLDYLDP ACFMHSLH
Source	in vitro E.coli expression system
Target Names	PRLR
Protein Names	Recommended name: Prolactin receptor Short name= PRL-R
Expression Region	25-616
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.